

Antibodies Product Datasheet

Product Name:	Nibrin Antibody
Catalog Number:	R12-3092
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	Nibrin Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Nibrin.
Concentration:	1 mg/ml
Gene Name:	NBN
Other Name:	Cell cycle regulatory protein P95; NBN; NBS; NIBRIN; NIJMEGEN BREAKAGE syndrome protein 1; p95-NBS1
Human Swiss Prot No.:	O60934
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Nibrin Antibodyantibody detects endogenous levels of Nibrin.
Storage and Stability:	Stable at -20°C for at least 1 year.