

Antibodies Product Datasheet

Product Name:	REEP5 Antibody
Catalog Number:	C51170
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C-terminal REEP5
Concentration:	1 mg/ml
Gene Name:	REEP5 C5orf18 DP1 TB2
Protein Name:	Receptor expression-enhancing protein 5 (Polyposis locus protein 1) (Protein TB2)
Other Name:	Receptor expression-enhancing protein 5 (Polyposis locus protein 1) (Protein TB2)
Human Swiss Prot No.:	Q00765
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of REEP5 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year