

Antibodies Product Datasheet

Product Name:	BM28 Antibody
Catalog Number:	C0259
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MCM2. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	MCM2
Protein Name:	DNA replication licensing factor MCM2
Other Name:	MCM2; BM28; CCNL1; CDCL1; KIAA0030; DNA replication licensing factor MCM2; Minichromosome maintenance protein 2 homolog; Nuclear protein BM28
Human Swiss Prot No.:	P49736
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BM28 Polyclonal Antibody detects endogenous levels of BM28 protein.
Storage and Stability:	-20°C/1 year