

Antibodies Product Datasheet

Product Name:	UMPS Antibody
Catalog Number:	C48471
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human UMPS
Concentration:	1 mg/ml
Gene Name:	UMPS OK/SW-cl.21
Protein Name:	Uridine 5'-monophosphate synthase (UMP synthase) [Includes: Orotate phosphoribosyltransferase (OPRT) (OPRTase) (EC 2.4.2.10)Orotidine 5'-phosphate decarboxylase (ODC) (EC 4.1.1.23) (OMPdecase)]
Other Name:	
Human Swiss Prot No.:	P11172
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of UMPS at Human, Mouse
Storage and Stability:	Store at -20° C/1 year