

Antibodies Product Datasheet

Product Name:	Tyrosine Hydrolase Mouse mAb(4B1)
Catalog Number:	N1457
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Tyrosine Hydrolase
Concentration:	1 mg/ml
Gene Name:	TH TYH
Protein Name:	Tyrosine Hydrolase
Other Name:	Tyrosine 3-monooxygenase (EC 1.14.16.2) (Tyrosine 3-hydroxylase) (TH)
Human Swiss Prot No.:	P07101
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.30% sodium azide.
Specificity:	This antibody detects endogenous levels of Tyrosine Hydrolase at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year