

Antibodies Product Datasheet

Product Name:	COL13A1 Antibody
Catalog Number:	C12217
Reactivity:	Human, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Collagen XIII alpha1. AA range:641-690
Concentration:	1 mg/ml
Gene Name:	COL13A1
Protein Name:	Collagen alpha-1(XIII) chain
Other Name:	COL13A1; Collagen alpha-1(XIII) chain; COLXIII A1
Human Swiss Prot No.:	Q5TAT6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	COL13A1 Polyclonal Antibody detects endogenous levels of COL13A1 protein.
Storage and Stability:	-20°C/1 year