

Antibodies Product Datasheet

Product Name:	mGluR-2/3 Antibody
Catalog Number:	C0207
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human mGluR2/3. AA range:823-872
Concentration:	1 mg/ml
Gene Name:	GRM2
Protein Name:	Metabotropic glutamate receptor 2
Other Name:	GRM2; GPRC1B; MGLUR2; Metabotropic glutamate receptor 2; mGluR2; GRM3; GPRC1C; MGLUR3; Metabotropic glutamate receptor 3; mGluR3
Human Swiss Prot No.:	Q14416/Q14832
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	mGluR-2/3 Polyclonal Antibody detects endogenous levels of mGluR-2/3 protein.
Storage and Stability:	-20°C/1 year