

Antibodies Product Datasheet

Product Name:	Hexokinase II mouse mAb
Catalog Number:	N0864
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Hexokinase II protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	hk2
Protein Name:	
Other Name:	DKFZp686M1669 ;Hexokinase 2 ;Hexokinase 2 muscle ;Hexokinase type II ;Hexokinase-2 ;HK 2 ;HK II ;HK2 ;HKII ;HxK 2 ;HxK2 ;HXK2_HUMAN ;Muscle form hexokinase .
Human Swiss Prot No.:	P52789
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Hexokinase II and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year