

Antibodies Product Datasheet

Product Name:	MMP-10 Antibody
Catalog Number:	C0263
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MMP-10. AA range:361-410
Concentration:	1 mg/ml
Gene Name:	MMP10
Protein Name:	Stromelysin-2
Other Name:	MMP10; STMY2; Stromelysin-2; SL-2; Matrix metalloproteinase-10; MMP-10; Transin-2
Human Swiss Prot No.:	P09238
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MMP-10 Polyclonal Antibody detects endogenous levels of MMP-10 protein.
Storage and Stability:	-20°C/1 year