

Antibodies Product Datasheet

Product Name:	BAF53 Antibody
Catalog Number:	C10255
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ACTL6A. AA range:201-250
Concentration:	1 mg/ml
Gene Name:	ACTL6A
Protein Name:	Actin-like protein 6A
Other Name:	ACTL6A; BAF53; BAF53A; INO80K; Actin-like protein 6A; 53 kDa BRG1-associated factor A; Actin-related protein Baf53a; ArpNbeta; BRG1-associated factor 53A; BAF53A; INO80 complex subunit K
Human Swiss Prot No.:	O96019
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BAF53 Polyclonal Antibody detects endogenous levels of BAF53 protein.
Storage and Stability:	-20°C/1 year