

Antibodies Product Datasheet

Product Name:	Caspase-9 mouse mAb
Catalog Number:	N0962
Reactivity:	Human, Transfected
Applications:	WB
Dilutions:	wb dilution 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Caspase-9 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	cas9
Protein Name:	Caspase9
Other Name:	APAF 3; APAF-3; APAF3; Apoptosis related cysteine peptidase; Apoptotic protease activating factor 3; Apoptotic protease Mch 6; Apoptotic protease Mch-6; Apoptotic protease-activating factor 3; CASP 9; CASP-9; CASP9; CASP9_HUMAN; Caspase 9 apoptosis related cysteine p
Human Swiss Prot No.:	P55211
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Caspase-9 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year