

## Antibodies Product Datasheet

<b>Product Name:</b>	ATF7IP Antibody
<b>Catalog Number:</b>	C49770
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal, Rabbit,IgG
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human ATF7IP
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ATF7IP MCAF MCAF1
<b>Protein Name:</b>	Activating transcription factor 7-interacting protein 1 (ATF-interacting protein) (ATF-IP) (ATF7-interacting protein) (ATFa-associated modulator) (hAM) (MBD1-containing chromatin-associated factor 1) (P621)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q6VMQ6
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of ATF7IP at Human, Mouse
<b>Storage and Stability:</b>	Store at -20° C/1 year