

Antibodies Product Datasheet

Product Name:	PSG11 Antibody
Catalog Number:	C50961
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PSG11
Concentration:	1 mg/ml
Gene Name:	PSG11 PSG13 PSG14
Protein Name:	Pregnancy-specific beta-1-glycoprotein 11 (PS-beta-G-11) (PSBG-11) (Pregnancy-specific glycoprotein 11) (Pregnancy-specific beta-1-glycoprotein 13) (PS-beta-G-13) (PSBG-13) (Pregnancy-specific glycoprotein 13)
Other Name:	
Human Swiss Prot No.:	Q9UQ72
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PSG11 at Human
Storage and Stability:	Store at -20° C/1 year