

Antibodies Product Datasheet

Product Name:	I-TAC Antibody
Catalog Number:	C40953
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human I-TAC
Concentration:	1 mg/ml
Gene Name:	CXCL11 ITAC SCYB11 SCYB9B
Protein Name:	I-TAC
Other Name:	C-X-C motif chemokine 11 (Beta-R1) (H174) (Interferon gamma-inducible protein 9) (IP-9) (Interferon-inducible T-cell alpha chemoattractant) (I-TAC) (Small-inducible cytokine B11)
Human Swiss Prot No.:	O14625
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human I-TAC
Storage and Stability:	-20°C/1 year