

Antibodies Product Datasheet

Product Name:	NMNAT1 Antibody
Catalog Number:	C48952
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human NMNAT1
Concentration:	1 mg/ml
Gene Name:	NMNAT1 NMNAT
Protein Name:	Nicotinamide mononucleotide adenylyltransferase 1 (NMN adenylyltransferase 1) (EC 2.7.7.1) (Nicotinate-nucleotide adenylyltransferase 1) (NaMN adenylyltransferase 1) (EC 2.7.7.18)
Other Name:	
Human Swiss Prot No.:	Q9HAN9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of NMNAT1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year