

Antibodies Product Datasheet

Product Name:	PRPS1L1 Antibody
Catalog Number:	C50633
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PRPS1L1
Concentration:	1 mg/ml
Gene Name:	PRPS1L1 PRPS3 PRPSL
Protein Name:	Ribose-phosphate pyrophosphokinase 3 (EC 2.7.6.1) (Phosphoribosyl pyrophosphate synthase 1-like 1) (PRPS1-like 1) (Phosphoribosyl pyrophosphate synthase III) (PRS-III)
Other Name:	
Human Swiss Prot No.:	P21108
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PRPS1L1 at Human
Storage and Stability:	Store at -20° C/1 year