

## Antibodies Product Datasheet

<b>Product Name:</b>	ILK (Phospho-Ser246) Antibody
<b>Catalog Number:</b>	A8115
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human ILK around the phosphorylation site of Ser246. AA range:212-261
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ILK
<b>Protein Name:</b>	Integrin-linked protein kinase
<b>Other Name:</b>	ILK; ILK1; ILK2; Integrin-linked protein kinase; 59 kDa serine/threonine-protein kinase; ILK-1; ILK-2; p59ILK
<b>Human Swiss Prot No.:</b>	Q13418
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Phospho-ILK (S246) Polyclonal Antibody detects endogenous levels of ILK protein only when phosphorylated at S246.
<b>Storage and Stability:</b>	-20°C/1 year