

Antibodies Product Datasheet

Product Name:	Separase Antibody
Catalog Number:	B0086
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SEPARASE. AA range:767-816
Concentration:	1 mg/ml
Gene Name:	ESPL1
Protein Name:	Separin
Other Name:	ESPL1; ESP1; KIAA0165; Separin; Caspase-like protein ESPL1; Extra spindle poles-like 1 protein; Separase
Human Swiss Prot No.:	Q14674
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Separase Polyclonal Antibody detects endogenous levels of Separase protein.
Storage and Stability:	-20°C/1 year