

Antibodies Product Datasheet

Product Name:	RUNX1T1 Antibody
Catalog Number:	C48992
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human RUNX1T1
Concentration:	1 mg/ml
Gene Name:	RUNX1T1 AML1T1 CBFA2T1 CDR ETO MTG8 ZMYND2
Protein Name:	Protein CBFA2T1 (Cyclin-D-related protein) (Eight twenty one protein) (Protein ETO) (Protein MTG8) (Zinc finger MYND domain-containing protein 2)
Other Name:	
Human Swiss Prot No.:	Q06455
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RUNX1T1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year