

Antibodies Product Datasheet

Product Name:	Cdk7 Antibody
Catalog Number:	C0353
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDK7. AA range:291-340
Concentration:	1 mg/ml
Gene Name:	CDK7
Protein Name:	Cyclin-dependent kinase 7
Other Name:	CDK7; CAK; CAK1; CDKN7; MO15; STK1; Cyclin-dependent kinase 7; 39 kDa protein kinase; p39 Mo15; CDK-activating kinase 1; Cell division protein kinase 7; Serine/threonine-protein kinase 1; TFIIF basal transcription factor complex kinase subu
Human Swiss Prot No.:	P50613
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdk7 Polyclonal Antibody detects endogenous levels of Cdk7 protein.
Storage and Stability:	-20°C/1 year