

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | LMTK2 Antibody |
| Catalog Number: | R12-2973 |
| Reactivity: | Human |
| Applications: | WB, IHC, ELISA |
| Dilutions: | LMTK2 Antibody |
| Modification: | Unmodified/Total |
| Source: | Rabbit Polyclonal |
| Purification: | The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human LMTK2. |
| Concentration: | 1 mg/ml |
| Gene Name: | LMTK2 |
| Other Name: | Lemur tyrosine kinase 2; Serine/threonine protein kinase KPI-2; Kinase/phosphatase/inhibitor 2; Apoptosis-associated tyrosine kinase 2; Brain-enriched kinase; hBREK; CDK5/p35-regulated kinase; CPRK Q8IWU2 |
| Human Swiss Prot No.: | Q8IWU2 |
| Formulation: | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Specificity: | LMTK2 Antibody antibody detects endogenous levels of LMTK2. |
| Storage and Stability: | Stable at -20°C for at least 1 year. |