

Antibodies Product Datasheet

Product Name:	p53(C-term) mouse mAb
Catalog Number:	N0793
Reactivity:	Human, Mouse, Monkey
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human p53 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	tp53
Protein Name:	
Other Name:	Antigen NY-CO-13; BCC7; Cellular tumor antigen p53; Cys 51 Stop; FLJ92943; HGNC11998; LFS1; p53; p53 Cellular Tumor Antigen; p53 tumor suppressor; P53_HUMAN; Phosphoprotein p53; Tp53; Transformation related protein 53; TRP53; Tumor protein p53; Tumor suppressor p53; Tum
Human Swiss Prot No.:	P04637
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of p53 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year