Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: Caspase-3 (Phospho-Ser150) Antibody

Catalog Number: A0835

Reactivity: Human, Mouse, Rat

Applications: WB, IHC-p, IF(paraffin section), ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000. Not yet tested in other applications.

Modification: Phospho Specific

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human Caspase 3 around the phosphorylation site of Ser150. AA

range:116-165

Concentration: 1 mg/ml

Gene Name: CASP3

Protein Name: Caspase3

Other Name: CASP3; CPP32; Caspase-3; CASP-3; Apopain; Cysteine protease

CPP32; CPP-32; Protein Yama; SREBP cleavage activity 1; SCA-1

Human Swiss Prot No.: P42574

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: Phospho-Caspase-3 (S150) Polyclonal Antibody detects endogenous

levels of Caspase-3 protein only when phosphorylated at S150.

Storage and Stability: -20°C/1 year