

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | HELLS Antibody |
| Catalog Number: | C48414 |
| Reactivity: | Human, Mouse |
| Applications: | WB |
| Dilutions: | WB 1:500-2000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal, Rabbit, IgG |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Synthesized peptide derived from human HELLS |
| Concentration: | 1 mg/ml |
| Gene Name: | HELLS PASG SMARCA6 Nbla10143 |
| Protein Name: | Lymphoid-specific helicase (EC 3.6.4.-) (Proliferation-associated SNF2-like protein) (SWI/SNF2-related matrix-associated actin-dependent regulator of chromatin subfamily A member 6) |
| Other Name: | |
| Human Swiss Prot No.: | Q9NRZ9 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | This antibody detects endogenous levels of HELLS at Human, Mouse |
| Storage and Stability: | Store at -20° C/1 year |