

Antibodies Product Datasheet

Product Name:	TRPM5 Antibody
Catalog Number:	C48601
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human TRPM5
Concentration:	1 mg/ml
Gene Name:	TRPM5 LTRPC5 MTR1
Protein Name:	Transient receptor potential cation channel subfamily M member 5 (Long transient receptor potential channel 5) (LTrpC-5) (LTrpC5) (MLSN1- and TRP-related gene 1 protein)
Other Name:	
Human Swiss Prot No.:	Q9NZQ8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of TRPM5 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year