

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | Syntenin-1 Monoclonal Antibody |
| Catalog Number: | N0752 |
| Reactivity: | Human,Rat,Dog,Pig |
| Applications: | WB |
| Dilutions: | Western Blot: 1/1000 - 1/2000. Not yet tested in other applications. |
| Modification: | Unmodified/ Total |
| Source: | Monoclonal Mouse |
| Purification: | Affinity purification |
| Immunogen: | Purified recombinant human Syntenin-1 protein fragments expressed in E.coli. |
| Concentration: | 1 mg/ml |
| Gene Name: | SDCBP |
| Protein Name: | Syntenin-1 |
| Other Name: | SDCBP; MDA9; SYCL; Syntenin-1; Melanoma differentiation-associated protein 9; MDA-9; Pro-TGF-alpha cytoplasmic domain-interacting protein 18; TACIP18; Scaffold protein Pbp1; Syndecan-binding protein 1 |
| Human Swiss Prot No.: | O00560 |
| Formulation: | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol. |
| Specificity: | Syntenin-1 Monoclonal Antibody detects endogenous levels of Syntenin-1 protein. |
| Storage and Stability: | -20°C/1 year |