

Antibodies Product Datasheet

Product Name:	EFNA1 (Cleaved-Ser182) Antibody
Catalog Number:	L0398
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EFNA1 (Cleaved-Ser182)
Concentration:	1 mg/ml
Gene Name:	EFNA1 EPLG1 LERK1 TNFAIP4
Protein Name:	EFNA1 (Cleaved-Ser182)
Other Name:	Ephrin-A1 (EPH-related receptor tyrosine kinase ligand 1) (LERK-1) (Immediate early response protein B61) (Tumor necrosis factor alpha-induced protein 4) (TNF alpha-induced protein 4) [Cleaved into: Ephrin-A1, secreted form]
Human Swiss Prot No.:	P20827
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat EFNA1 (Cleaved-Ser182, protein was cleaved amino acid sequence between 182-183)
Storage and Stability:	-20°C/1 year