
Antibodies Product Datasheet

Product Name:	Cyclin B1 mouse Monoclonal Antibody(5C1)
Catalog Number:	N1643
Reactivity:	Human,Rat,Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Cyclin B1 at AA range of 60-140
Concentration:	1 mg/ml
Gene Name:	CCNB1
Protein Name:	CCNB1
Other Name:	CCNB1
Human Swiss Prot No.:	P14635
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cyclin B1 protein detects endogenous levels of CCNB1
Storage and Stability:	-20°C/1 year