

Antibodies Product Datasheet

Product Name:	RNF130 Antibody
Catalog Number:	C20560
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RNF130. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	RNF130
Protein Name:	E3 ubiquitin-protein ligase RNF130
Other Name:	RNF130; E3 ubiquitin-protein ligase RNF130; Goliath homolog; H-Goliath; RING finger protein 130
Human Swiss Prot No.:	Q86XS8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RNF130 Polyclonal Antibody detects endogenous levels of RNF130 protein.
Storage and Stability:	-20°C/1 year