

Antibodies Product Datasheet

Product Name:	PKC β (Phospho-Ser661) Antibody
Catalog Number:	A8169
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PKCB around the phosphorylation site of Ser661. AA range:622-671
Concentration:	1 mg/ml
Gene Name:	PRKCB
Protein Name:	Protein kinase C beta type
Other Name:	PRKCB; PKCB; PRKCB1; Protein kinase C beta type; PKC-B; PKC-beta
Human Swiss Prot No.:	P05771
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-PKC β (S661) Polyclonal Antibody detects endogenous levels of PKC β protein only when phosphorylated at S661.
Storage and Stability:	-20°C/1 year