

## **Assay Biotechnology Company**

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

## **Antibodies Product Datasheet**

**Product Name:** AMPKα1 Antibody

Catalog Number: B7007

Reactivity: Human, Mouse, Rat

**Applications:** WB, IHC-p, IF/ICC, ELISA

**Dilutions:** Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in

other applications.

**Modification:** Unmodified/Total

**Source:** Polyclonal Rabbit

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Immunogen:** The antiserum was produced against synthesized peptide derived from

human AMPK1. AA range:451-500

Concentration: 1 mg/ml

Gene Name: PRKAA1

**Protein Name:** 5'-AMP-activated protein kinase catalytic subunit alpha-1

Other Name: PRKAA1; AMPK1; 5'-AMP-activated protein kinase catalytic subunit

alpha-1; AMPK subunit alpha-1; Acetyl-CoA carboxylase kinase; ACACA kinase; Hydroxymethylglutaryl-CoA reductase kinase; HMGCR kinase;

Tau-protein kinase PRKAA1

Human Swiss Prot No.: Q13131

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: AMPKα1 Polyclonal Antibody detects endogenous levels of AMPKα1

protein.

Storage and Stability: -20°C/1 year