

Antibodies Product Datasheet

Product Name:	MAK (Phospho-Tyr159) Antibody
Catalog Number:	A8128
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAK around the phosphorylation site of Tyr159. AA range:126-175
Concentration:	1 mg/ml
Gene Name:	MAK
Protein Name:	Serine/threonine-protein kinase MAK
Other Name:	MAK; Serine/threonine-protein kinase MAK; Male germ cell-associated kinase
Human Swiss Prot No.:	P20794
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-MAK (Y159) Polyclonal Antibody detects endogenous levels of MAK protein only when phosphorylated at Y159.
Storage and Stability:	-20°C/1 year