

## Antibodies Product Datasheet

<b>Product Name:</b>	Survivin Monoclonal Antibody(Mix)
<b>Catalog Number:</b>	N1515
<b>Reactivity:</b>	Human,Rat
<b>Applications:</b>	WB,IHC-p,IF(paraffin section)
<b>Dilutions:</b>	WB: 1:1000-2000 IHC: 1:200-500
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Recombinant Protein of Survivin
<b>Concentration:</b>	
<b>Gene Name:</b>	BIRC5
<b>Protein Name:</b>	Baculoviral IAP repeat-containing protein 5
<b>Other Name:</b>	BIRC5; API4; IAP4; Baculoviral IAP repeat-containing protein 5; Apoptosis inhibitor 4; Apoptosis inhibitor survivin
<b>Human Swiss Prot No.:</b>	O15392
<b>Formulation:</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous Survivin protein.
<b>Storage and Stability:</b>	-20°C/1 year