

Antibodies Product Datasheet

Product Name:	Nek9 Antibody
Catalog Number:	B1159
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NEK9. AA range:176-225
Concentration:	1 mg/ml
Gene Name:	NEK9
Protein Name:	Serine/threonine-protein kinase Nek9
Other Name:	NEK9; KIAA1995; NEK8; NERCC; Serine/threonine-protein kinase Nek9; Nercc1 kinase; Never in mitosis A-related kinase 9; NimA-related protein kinase 9; NimA-related kinase 8; Nek8
Human Swiss Prot No.:	Q8TD19
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Nek9 Polyclonal Antibody detects endogenous levels of Nek9 protein.
Storage and Stability:	-20°C/1 year