
Antibodies Product Datasheet

Product Name:	EGF Receptor (Phospho-Ser1036/1037) Antibody
Catalog Number:	A9495
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human EGF Receptor (Ser1036 and 1037)
Concentration:	1 mg/ml
Gene Name:	EGFR ERBB ERBB1 HER1
Protein Name:	EGF Receptor (Ser1036/1037)
Other Name:	Epidermal growth factor receptor (EC 2.7.10.1) (Proto-oncogene c-ErbB-1) (Receptor tyrosine-protein kinase erbB-1)
Human Swiss Prot No.:	P00533
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat EGF Receptor (phospho-Ser1036 or 1037)
Storage and Stability:	-20°C/1 year