

Antibodies Product Datasheet

Product Name:	Endothelin B Receptor Antibody
Catalog Number:	C40189
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Endothelin B Receptor
Concentration:	1 mg/ml
Gene Name:	EDNRB
Protein Name:	Endothelin B receptor (ET-B) (ET-BR) (Endothelin receptor non-selective type)
Other Name:	Endothelin B receptor (ET-B) (ET-BR) (Endothelin receptor non-selective type); YM3547
Human Swiss Prot No.:	P24530
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Endothelin B Receptor protein(A221) detects endogenous levels of Endothelin B Receptor
Storage and Stability:	-20°C/1 year