
Antibodies Product Datasheet

Product Name:	Insulin R Antibody
Catalog Number:	B0492
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IR. AA range:881-930
Concentration:	1 mg/ml
Gene Name:	INSR
Protein Name:	Insulin receptor
Other Name:	INSR; Insulin receptor; IR; CD antigen CD220
Human Swiss Prot No.:	P06213
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Insulin R Polyclonal Antibody detects endogenous levels of Insulin R protein.
Storage and Stability:	-20°C/1 year