

Antibodies Product Datasheet

Product Name:	PSMB9 Antibody
Catalog Number:	C48271
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PSMB9
Concentration:	1 mg/ml
Gene Name:	PSMB9 LMP2 PSMB6i RING12
Protein Name:	Proteasome subunit beta type-9 (EC 3.4.25.1) (Low molecular mass protein 2) (Macropain chain 7) (Multicatalytic endopeptidase complex chain 7) (Proteasome chain 7) (Proteasome subunit beta-1i) (Really interesting new gene 12 protein)
Other Name:	
Human Swiss Prot No.:	P28065
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PSMB9 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year