

## Antibodies Product Datasheet

<b>Product Name:</b>	CK19 Monoclonal Antibody(11F5)
<b>Catalog Number:</b>	N1315
<b>Reactivity:</b>	Human,Mouse,Rat
<b>Applications:</b>	IHC-p,IF(paraffin section)
<b>Dilutions:</b>	IHC: 1:200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthetic Peptide of CK19
<b>Concentration:</b>	
<b>Gene Name:</b>	KRT19
<b>Protein Name:</b>	Keratin type I cytoskeletal 19
<b>Other Name:</b>	KRT19; Keratin, type I cytoskeletal 19; Cytokeratin-19; CK-19; Keratin-19; K19
<b>Human Swiss Prot No.:</b>	P08727
<b>Formulation:</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous CK19 proteins.
<b>Storage and Stability:</b>	-20°C/1 year