

Antibodies Product Datasheet

Product Name:	Melan-A Antibody
Catalog Number:	C0257
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MART-1. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	MLANA
Protein Name:	Melanoma antigen recognized by T-cells 1
Other Name:	MLANA; MART1; Melanoma antigen recognized by T-cells 1; MART-1; Antigen LB39-AA; Antigen SK29-AA; Protein Melan-A
Human Swiss Prot No.:	Q16655
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Melan-A Polyclonal Antibody detects endogenous levels of Melan-A protein.
Storage and Stability:	-20°C/1 year