

Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: FoxO3A Antibody

Catalog Number: B7087

Reactivity: Human, Mouse, Rat

Applications: WB, IHC-p, IF/ICC, ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in

other applications.

Modification: Unmodified/Total

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human FKHRL1. AA range:220-269

Concentration: 1 mg/ml

Gene Name: FOXO3

Protein Name: Forkhead box protein O3

Other Name: FOXO3; FKHRL1; FOXO3A; Forkhead box protein O3; AF6q21 protein;

Forkhead in rhabdomyosarcoma-like 1

Human Swiss Prot No.: 043524

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: FoxO3A Polyclonal Antibody detects endogenous levels of FoxO3A

protein.

Storage and Stability: -20°C/1 year