

Antibodies Product Datasheet

Product Name:	COP1 Antibody
Catalog Number:	C15565
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RFWD2. AA range:661-710
Concentration:	1 mg/ml
Gene Name:	RFWD2
Protein Name:	E3 ubiquitin-protein ligase RFWD2
Other Name:	RFWD2; COP1; RNF200; E3 ubiquitin-protein ligase RFWD2; Constitutive photomorphogenesis protein 1 homolog; hCOP1; RING finger and WD repeat domain protein 2; RING finger protein 200
Human Swiss Prot No.:	Q8NHY2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	COP1 Polyclonal Antibody detects endogenous levels of COP1 protein.
Storage and Stability:	-20°C/1 year