
Antibodies Product Datasheet

Product Name:	CYP1A2 Antibody
Catalog Number:	C12250
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 1A2. AA range:331-380
Concentration:	1 mg/ml
Gene Name:	CYP1A2
Protein Name:	Cytochrome P450 1A2
Other Name:	CYP1A2; Cytochrome P450 1A2; CYP1A2; Cytochrome P(3)450; Cytochrome P450 4; Cytochrome P450-P3
Human Swiss Prot No.:	P05177
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYP1A2 Polyclonal Antibody detects endogenous levels of CYP1A2 protein.
Storage and Stability:	-20°C/1 year