

Antibodies Product Datasheet

Product Name:	GGCX Antibody
Catalog Number:	C48406
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GGCX
Concentration:	1 mg/ml
Gene Name:	GGCX GC
Protein Name:	Vitamin K-dependent gamma-carboxylase (EC 4.1.1.90) (Gamma-glutamyl carboxylase) (Peptidyl-glutamate 4-carboxylase) (Vitamin K gamma glutamyl carboxylase)
Other Name:	
Human Swiss Prot No.:	P38435
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GGCX at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year