

Antibodies Product Datasheet

Product Name:	PKC ϵ (Phospho-Ser729) Antibody
Catalog Number:	A0802
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PKC epsilon around the phosphorylation site of Ser729. AA range:688-737
Concentration:	1 mg/ml
Gene Name:	PRKCE
Protein Name:	Protein kinase C epsilon type
Other Name:	PRKCE; PKCE; Protein kinase C epsilon type; nPKC-epsilon
Human Swiss Prot No.:	Q02156
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-PKC ϵ (S729) Polyclonal Antibody detects endogenous levels of PKC ϵ protein only when phosphorylated at S729.
Storage and Stability:	-20°C/1 year