
Antibodies Product Datasheet

Product Name:	β -Tubulin Monoclonal Antibody(5G3)
Catalog Number:	N1294
Reactivity:	Human,Rat,Mouse,Monkey,Dog,Chicken,Hamster,Rabbit,Sheep,Insect,Yeast,Fish
Applications:	WB,IF/ICC,IHC-p
Dilutions:	WB: 1:5000 IHC: 1:200 IF 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of β -Tubulin
Concentration:	1 mg/ml
Gene Name:	TUBB3
Protein Name:	Tubulin beta-3 chain
Other Name:	TUBB3; TUBB4; Tubulin beta-3 chain; Tubulin beta-4 chain; Tubulin beta-III
Human Swiss Prot No.:	Q13509
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous β -tubulin protein.
Storage and Stability:	-20°C/1 year