

Antibodies Product Datasheet

Product Name:	AurB/C Antibody
Catalog Number:	R12-2525
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	AurB/C Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AurB/C.
Concentration:	1 mg/ml
Gene Name:	AURKB/AURKC
Other Name:	AIK2; AIM-1; AIM1; ARK2; AURKB; Aurora- and Ipl1-like midbody-associated protein 1; Aurora-B; Aurora-related kinase 2; Aurora/IPL1-related kinase 2; Serine/threonine protein kinase 12; STK-1; STK12; AIE1; AIE2; AURKC; aurora-C; aurora/Ipl1-related kinase
Human Swiss Prot No.:	Q96GD4/Q9UQB9
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	AurB/C Antibody antibody detects endogenous levels of AurB/C.
Storage and Stability:	Stable at -20°C for at least 1 year.