

Antibodies Product Datasheet

Product Name:	Caspase-6 (Phospho-Ser257) Antibody
Catalog Number:	A0058
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Caspase 6 around the phosphorylation site of Ser257. AA range:223-272
Concentration:	1 mg/ml
Gene Name:	CASP6
Protein Name:	Caspase6
Other Name:	CASP6; MCH2; Caspase-6; CASP-6; Apoptotic protease Mch-2
Human Swiss Prot No.:	P55212
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Caspase-6 (S257) Polyclonal Antibody detects endogenous levels of Caspase-6 protein only when phosphorylated at S257.
Storage and Stability:	-20°C/1 year