

Antibodies Product Datasheet

Product Name:	Cyclin D1 Antibody
Catalog Number:	B0417
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cyclin D1. AA range:246-295
Concentration:	1 mg/ml
Gene Name:	CCND1
Protein Name:	G1/S-specific cyclin-D1
Other Name:	CCND1; BCL1; PRAD1; G1/S-specific cyclin-D1; B-cell lymphoma 1
Human Swiss Prot No.:	protein; BCL-1; BCL-1 oncogene; PRAD1 oncogene P24385
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cyclin D1 Polyclonal Antibody detects endogenous levels of Cyclin D1 protein.
Storage and Stability:	-20°C/1 year