

Antibodies Product Datasheet

Product Name:	GDPD2 Antibody
Catalog Number:	C50382
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GDPD2
Concentration:	1 mg/ml
Gene Name:	GDPD2 GDE3 OBDPF UNQ1935/PRO4418
Protein Name:	Glycerophosphoinositol inositolphosphodiesterase GDPD2 (EC 3.1.4.43) (Glycerophosphodiester phosphodiesterase 3) (Glycerophosphodiester phosphodiesterase domain-containing protein 2) (Osteoblast differentiation promoting factor)
Other Name:	
Human Swiss Prot No.:	Q9HCC8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GDPD2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year